ISO 11140-4
Compliant

Tested by 3rd
Party Test Institute

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
User-friendly & Comprehensive Software Package

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Utilize the Full Potential of SteriSense

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
Operation:	Autoclave
Process:	Steam Sterilization
Product:	Axi
Version:	Autoclave 1
Session Start:	26-03-2019 11:49:20
Session Stop:	26-03-2019 12:41:34
Operator:	autoclave
Batch Number:	1
Test Result:	Passed

ISO 11140-4/17665 Electronic Bowie & Dick	
Item:	ISO 11140-4/17665 Electronic Bowie & Dick Test 1
Category:	Chemical
Batch Number:	1
Test standard status:	ISO 17665 compliant
Total Test Result:	Passed

Input parameters	
Process Start:	26-03-2019 11:49:20
Process Stop:	26-03-2019 12:41:34
Operator:	autoclave
Improvement Type:	BI-Pho
Table name:	Report: BI-Pho, Autoclave

Process plan	
Sterilization validation method:	Original sample culture
Sterilization temperature setpoint:	134.00 °C
Upper Offset:	0.00 °C
Maximum temperature difference:	0.00 °C
Maximum temperature difference limit (°C):	5.00 °C
Sterilization duration:	00:03:00

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Technical Specification of SteriSense®

Up to 1,000 test
cycles between
each calibration

SteriSense Unit	
ISO Standard Compliance:	11140-4:2007 type B1, B2 and B3 at 134 °C
3rd Party Test Institute:	SAL GmbH
Type of Sterilizers:	Steam Sterilizers / Autoclaves qualified for sterilization at +134 °C according to EN 285 (volume larger than 60 L) and EN ISO 17665 (moist heat autoclaves)
Measuring Principle:	Piezoresistive / Electrical Resistance
Sensor Type:	Strain Gauge / Pt1000
Temperature Measuring Range:	0 to +140 °C (Calibrated +25 to +140 °C)
Temperature Accuracy: +25 to +140 °C:	± 0.05 °C
Pressure Measuring Range:	10 mBar to 6 Bar
Accuracy:	± 0.25% Full Scale (± 15 mBar)
Logger Housing Material:	316L Stainless Steel
PCD Housing Material:	PEEK
Operating Temperature:	-20 to +150 °C
Operating Pressure:	0.001 mBar to 10 Bar ABS
Diameter:	25 mm
Length:	125 mm including PCD (30 mm)
Weight with Battery:	240 g
Memory Capacity:	40,000 Data Points / 10,000 Samples
Minimum Sample Rate:	1 Second
Maximum Sample Rate:	24 Hours
Intrinsically Safe:	Ex II1GD ia IIC T3 Ga, -55 °C ≤ Tamb ≤ +105 °C
Time Accuracy:	± 5 Seconds per 24 Hours
Battery:	TSP 150L Battery
Expected Battery Lifetime:	1,000 Tests / 500 Hours (at +134 °C with Sample Rate of 1 Sec)
If equipment is used in ATEX environment, special conditions for safe use are stated in ATEX certificates , section 17 must be considered.	